



Vitrectomy Cutters - 25Ga

Vitrectomy Cutters, 25Ga

- 5157** 30 PSI **Alcon:** Accurus - 2500 cpm
DORC: Associate - 2500 cpm
- B+L:** Millennium, Premiere - 750 cpm
DORC: Harmony Total TTC - 1500 cpm
 5022 - Pressure Line Adapter required
- 5140** 25 PSI **Carl Zeiss Meditec/IOL Tech:** Pentasys - 1500 cpm
Carl Zeiss Meditec: Visalis 500
Fritz Ruck: Pentasys, Pentasys 2 - 1500 cpm
Croma Pharma/Corneal: Open Phaco - 1500 cpm
Optikon: Pulsar², Assistant - 2000 cpm
AMO: Gemini - 2000 cpm
Nidek: CV-24000 - 2500 cpm
Nidek: CV-12000, VT-5000 - 800 cpm



5022 - Pressure Line Adapter - 5 per pack
 Connects 5140 with DORC Harmony Total TTC



Fiber Optic Probes, 25Ga

- 5135** Fiber Optic Probe, 25Ga - Alcon Accurus
5258 Fiber Optic Probe, 25Ga - Alcon Constellation
5255 Wide Angle Fiber Optic Probe, 25Ga - Alcon Constellation
5138 Fiber Optic Probe, 25Ga - B+L Millennium
5259 Fiber Optic Probe, 25Ga - B+L Stellaris PC
5238 Wide Angle Fiber Optic Probe, 25Ga - B+L Stellaris PC
5126 Fiber Optic Probe, 25Ga - Ocutome connector
5249 Chandelier Fiber with Adjustable Ring, 25Ga/0.5mm - B+L Stellaris PC
5250 Chandelier Fiber with Adjustable Ring, 25Ga/0.5mm - Alcon Constellation
5228 Chandelier Fiber with Adjustable Ring, 25Ga/0.5mm - Ocutome connector




Silicone Tip Cannulas, 25Ga


- 5160** Brush Tip Cannula, 25Ga - 0.6mm silicone brush tip
5152 Soft Tip Cannula, 25Ga - 2.0mm silicone tip
5156 Charles Flute Cannula, 25Ga - Blunt tip



Brush Tip



Soft Tip



Charles Flute

One-Step Trocar System - 25Ga

Available January 2012

Second Generation

- Self-sealing wound construction
- Ultra sharp knife for smooth incision
- Unique closure valve system
- 3.0 + 4.0 mm scleral marker

5244 One-Step Trocar System, 25Ga


Each package includes:

- Three Ultra sharp Trocar Cannula Knives
- Three Preloaded self-sealing Trocar Cannulas
- One Universal Infusion Line




Trocar Cannula Knife

Knife Inserter consist of ultra sharp blade with preloaded Self-Sealing Trocar Cannula



Self-Sealing Trocar Cannula

Stainless steel Trocar Cannula is complete with self-sealing silicone valve, eliminating need for plugs. Valve offers superior sealing and smooth instruments extraction



Universal Infusion Line

- Instantly connects to Trocar Cannula
- Designed to offer maximum infusion flow



One-Step Trocar System - 25Ga

First Generation

- Ultra sharp blade mounted to plastic handle
- Low pressure insertion
- Preloaded Trocar cannula
- Self-sealing silicone valve

5150 One Step Trocar System, 25Ga

Each package includes:


- Four Trocar Cannula Inserters
- Four Preloaded Self-sealing Trocar Cannulas

Watch video clip on [MicroVision website](#)



Trocar Cannula Inserter

Trocar Cannula Inserter comes with Cannula pre-installed. Delivered in retracted position.



Self-Sealing Trocar Cannula

Stainless steel Trocar Cannula comes with self-sealing valve, eliminating need for plugs.

