


Iris Retractors

CE



- **Round polypropylene fiber**
- Ideal for temporary iris dilation
- Insertions through 0.5 mm incision
- Unique, easy to adjust flat silicone ring
- Floppy Iris Syndrome management in Flomax[®] patients



Disposable Iris Retractors

5016 -1	1 pack of 4 pcs - 4 retractors
5016 -5	5 packs of 4 pcs - 20 retractors
5067 - 1	1 pack of 5 pcs - 5 retractors
5067- 6	6 packs of 5 pcs - 30 retractors



Reusable Iris Retractors


5127 | 1 pack of 6 pcs - 6 retractors

Steam autoclavable container

Watch video clip on MicroVision website

Flatar[®] Iris Retractor

CE



- **Flat polypropylene fiber**
- Full control in placement and removal with flat fiber
- Larger Iris contact area - less chance of iris tissue tearing
- Less traumatic for iris tissue
- More stable than round Iris Retractor
- Unique, easy to adjust flat silicone ring
- Insertions through 0.5 mm incision



5213 | Flatar[®] Iris Retractor - 1 pack of 5 pcs