

# Posterior Vitrectomy Cutters, 20Ga

## Vitrectomy Cutters - 800 cuts per minute



### Disposable Cutters



**5064** Alcon: Accurus  
30 PSI  
DORC: Harmony Total TTC

**5020** Alcon: STTO  
20 PSI  
Nidek: CV-12000

**5099** AMO: Gemini  
25 PSI  
B&L: Premiere, Millennium, Daisy  
Optikon: Pulsar<sup>2</sup>, Assistant  
Carl Zeiss Meditec/IOL Tech: Pentasys  
Fritz Ruck: Pentasys  
Croma Pharma/Corneal: Open Phaco  
Nidek: CV-24000, VT-5000  
MID Labs: MVS, SupraVit  
Syntec: VitMan

### Reusable Cutters



**5065** Alcon: Accurus  
30 PSI  
DORC: Harmony Total TTC

**5026** Alcon: STTO  
20 PSI  
Nidek: CV-12000

**5100** AMO: Gemini  
25 PSI  
B&L: Premiere, Millennium, Daisy  
Optikon: Pulsar<sup>2</sup>, Assistant  
Carl Zeiss Meditec/IOL Tech: Pentasys  
Fritz Ruck: Pentasys  
Croma Pharma/Corneal: Open Phaco  
Nidek: CV-24000, VT-5000  
MID Labs: MVS, SupraVit  
Syntec: VitMan

## High Speed Vitrectomy Cutters

**5124** Alcon: Accurus - 2500 cpm  
30 PSI  
DORC: Associate - 2500 cpm

**5215** DORC: Harmony Total TTC - 1500 cpm  
25 PSI

**5168** Carl Zeiss Meditec/IOL Tech: Pentasys - 1500 cpm  
25 PSI  
Fritz Ruck: Pentasys - 1500 cpm  
Croma Pharma/Corneal: Open Phaco - 1500 cpm  
Optikon: Pulsar<sup>2</sup>, Assistant - 2000 cpm  
AMO: Gemini - 2000 cpm  
Nidek: CV-24000 - 2500 cpm

