

# Anterior Vitrectomy Cutters - 20Ga

## Vitrectomy Cutters - 800 cuts per minute

### Disposable Cutters with Irrigation



**5045**  
20 PSI

**Alcon:** Legacy, Universal  
**AMO:** Sovereign, Signature, AMO Plus, Prestige  
**Staar Surgical:** Wave  
**Surgical Design:** Ocusystem  
**American Optisurgical:** Horizon, Vizual  
**Nidek:** CV-12000, CV-6000, CV-7000

**5118**  
25 PSI

**AMO:** Gemini  
**B+L:** Stellaris, Millennium, Premiere, Protégé  
**Optikon:** Pulsar<sup>2</sup>, Assistant  
**Croma Pharma/Corneal:** Open Phaco  
**Carl Zeiss Meditec/IOL Tech/Fritz Ruck:** Pentasys  
**Carl Zeiss:** Visalis 100  
**Medical Technical Products**  
**Mediphacos**  
**Nidek:** CV-24000

### Reusable Cutters with Irrigation



**5050**  
20 PSI

**Alcon:** Legacy, Universal  
**AMO:** Sovereign, Signature, AMO Plus, Prestige  
**Staar Surgical:** Wave  
**Surgical Design:** Ocusystem  
**American Optisurgical:** Horizon, Vizual  
**Nidek:** CV-12000, CV-6000, CV-7000

**5119**  
25 PSI

**AMO:** Gemini  
**B+L:** Stellaris, Millennium, Premiere, Protégé  
**Optikon:** Pulsar<sup>2</sup>, Assistant  
**Croma Pharma/Corneal:** Open Phaco  
**Carl Zeiss Meditec/IOL Tech/Fritz Ruck:** Pentasys  
**Carl Zeiss:** Visalis 100  
**Medical Technical Products**  
**Mediphacos**  
**Nidek:** CV-24000

### Disposable Cutters without Irrigation



**5068**  
20 PSI

**Alcon:** Legacy, Universal  
**AMO:** Sovereign, Signature, AMO Plus, Prestige  
**Staar Surgical:** Wave  
**Surgical Design:** Ocusystem  
**American Optisurgical:** Horizon, Vizual  
**Nidek:** CV-12000, CV-6000, CV-7000


**5101**  
25 PSI

**AMO:** Gemini  
**B+L:** Stellaris, Millennium, Premiere, Protégé  
**Optikon:** Pulsar<sup>2</sup>, Assistant  
**Croma Pharma/Corneal:** Open Phaco  
**Carl Zeiss Meditec/IOL Tech/Fritz Ruck:** Pentasys  
**Carl Zeiss:** Visalis 100  
**Medical Technical Products**  
**Mediphacos**  
**Nidek:** CV-24000

**5142**  
30 PSI

**Alcon:** Accurus, Infiniti  
**DORC:** Harmony Total TTC  
 Packaged with 21Ga Irrigation Cannula

### Reusable Cutters without Irrigation



**5070**  
20 PSI

**Alcon:** Legacy, Universal  
**AMO:** Sovereign, Signature, AMO Plus, Prestige  
**Staar Surgical:** Wave  
**Surgical Design:** Ocusystem  
**American Optisurgical:** Horizon, Vizual  
**Nidek:** CV-12000, CV-6000, CV-7000

**5102**  
25 PSI

**AMO:** Gemini  
**B+L:** Stellaris, Millennium, Premiere, Protégé  
**Optikon:** Pulsar<sup>2</sup>, Assistant  
**Croma Pharma/Corneal:** Open Phaco  
**Carl Zeiss Meditec/IOL Tech/Fritz Ruck:** Pentasys  
**Carl Zeiss:** Visalis 100  
**Medical Technical Products**  
**Mediphacos**  
**Nidek:** CV-24000

**5143**  
30 PSI

**Alcon:** Accurus, Infiniti  
**DORC:** Harmony Total TTC  
 Packaged with 21Ga Irrigation Cannula

### Disposable Irrigation Sleeve

**5035** - 5 per package



### Reusable Irrigation Sleeve

**5034** - 1 per package

