

Anterior Vitrectomy Cutters, 20Ga

Vitrectomy Cutters - 800 cuts per minute



Cutters with Irrigation



5045
20 PSI

Alcon: Legacy, Universal
AMO: Sovereign, Signature, AMO Plus, Prestige
Staar Surgical: Wave
Surgical Design: Ocusystem
American Optisurgical: Horizon, Vizual
Nidek: CV-12000, CV-6000, CV-7000

5118
25 PSI

AMO: Gemini
B&L: Stellaris, Millennium, Premiere, Protégé
Optikon: Pulsar², Assistant
Croma Pharma/Corneal: Open Phaco
Carl Zeiss Meditec/IOL Tech/Fritz Ruck: Pentasys
Carl Zeiss: Visalis 100
Medical Technical Products
Mediphacos
Nidek: CV-24000

Reusable Cutters with Irrigation



5050
20 PSI

Alcon: Legacy, Universal
AMO: Sovereign, Signature, AMO Plus, Prestige
Staar Surgical: Wave
Surgical Design: Ocusystem
American Optisurgical: Horizon, Vizual
Nidek: CV-12000, CV-6000, CV-7000

5119
25 PSI

AMO: Gemini
B&L: Stellaris, Millennium, Premiere, Protégé
Optikon: Pulsar², Assistant
Croma Pharma/Corneal: Open Phaco
Carl Zeiss Meditec/IOL Tech/Fritz Ruck: Pentasys
Carl Zeiss: Visalis 100
Medical Technical Products
Mediphacos
Nidek: CV-24000

Cutters without Irrigation



5068
20 PSI

Alcon: Legacy, Universal
AMO: Sovereign, Signature, AMO Plus, Prestige
Staar Surgical: Wave
Surgical Design: Ocusystem
American Optisurgical: Horizon, Vizual
Nidek: CV-12000, CV-6000, CV-7000


5101
25 PSI

AMO: Gemini
B&L: Stellaris, Millennium, Premiere, Protégé
Optikon: Pulsar², Assistant
Croma Pharma/Corneal: Open Phaco
Carl Zeiss Meditec/IOL Tech/Fritz Ruck: Pentasys
Carl Zeiss: Visalis 100
Medical Technical Products
Mediphacos
Nidek: CV-24000

5142
30 PSI

Alcon: Accurus, Infiniti
DORC: Harmony Total TTC
 Packaged with 21Ga Irrigation Cannula

Reusable Cutters without Irrigation



5070
20 PSI

Alcon: Legacy, Universal
AMO: Sovereign, Signature, AMO Plus, Prestige
Staar Surgical: Wave
Surgical Design: Ocusystem
American Optisurgical: Horizon, Vizual
Nidek: CV-12000, CV-6000, CV-7000

5102
25 PSI

AMO: Gemini
B&L: Stellaris, Millennium, Premiere, Protégé
Optikon: Pulsar², Assistant
Croma Pharma/Corneal: Open Phaco
Carl Zeiss Meditec/IOL Tech/Fritz Ruck: Pentasys
Carl Zeiss: Visalis 100
Medical Technical Products
Mediphacos
Nidek: CV-24000

5143
30 PSI

Alcon: Accurus, Infiniti
DORC: Harmony Total TTC
 Packaged with 21Ga Irrigation Cannula

Disposable Irrigation Sleeve

5035 - 5 per package



Reusable Irrigation Sleeve

5034 - 1 per package

